

**To include on grant or paper submissions:**

Cell lines were validated by STR DNA fingerprinting using the Promega 16 High Sensitivity STR Kit (Catalog # DC2100). The STR profiles were compared to online search databases (DSMZ/ATCC/JCRB/RIKEN) of approximately 2500 known profiles; along with the MD Anderson Characterized Cell Line Core (CCLC) database of approximately 2600 known profiles. The STR profiles matched known DNA fingerprints or were unique.

**In Acknowledgements:**

STR DNA fingerprinting was done by the Cytogenetics and Cell Authentication Core (CCAC) at M.D. Anderson Cancer Center.